



Suction Instruments



29-816
Universal cannula, double curved
230 mm



▲
Ø 7 mm
Chart. 21

POOL SUCTION TUBE

29-860
225 mm



▲
Ø 10 mm
Chart. 30

POOLE SUCTION TUBE

29-862 - (SMALL)
29-863 - (MEDIUM)
29-864 - (LARGE)



YANKAUER


29-944
with central bore, 2.0 mm diam.
285 mm



DE BAKEY SUCTION TUBE

29-946

with central bore, 3.0 mm diam. with lateral slots
275 mm



YANKAUER SUCTION TUBE

29-671

with especially large central bore, 4.8 mm diam.
285 mm



DE BAKEY SUCTION TUBE

29-948

with finger cut-off and stylet 270 mm



DE BAKEY SUCTION TUBE

29-947

with finger cut-off and stylet 160 mm



ADSON SUCTION TUBE

29-917
with finger cut-off
190 mm



ADSON SUCTION TUBE

29-918
with finger cut-off and stylet
165 mm



PLESTER SUTION TUBE

195 mm

	diam.
29-939	1,5 mm
29-940	2,0 mm
29-941	2,5 mm
29-942	3,0 mm
29-355	5,0 mm



FERGUSON SUCTION TUBE

90° to curve, 180 mm

	diam.
29-350	1,5 mm
29-351	2,0 mm
29-352	2,5 mm
29-353	3,0 mm
29-354	4,0 mm



olive diam. 6 - 9 mm

FERGUSSON SUCTION TUBE

29-933 Ø diam. 4,0 mm
90° to curve, 180 mm



olive diam. 6 - 9 mm

FRAZIER SUCTION TUBE

diam.
29-923 3,0 mm
29-925 5,0 mm
90° to curve, 180 mm



1/2

FRAZIER SUCTION TUBE

diam.
29-124 2,3 mm
29-125 2,7 mm
29-127 3,4 mm
29-129 4,0 mm
90° to curve, 110 mm



LEMPERT SUCTION TUBE

diam.
29-937 2,0 mm
29-938 3,0 mm
with finger cut-off and stylet
200 mm



FORMBY SUCTION TUBE

diam.
29-932 3,3 mm
210 mm



BARRON SUCTION TUBE

diam.
29-935 1,3 mm
29-936 2,0 mm
with finger cut-off and stylet
160 mm



SCHUHKNECHT SUCTION TUBE

diam. diam.
29-756 20 swg 1,0 mm
29-757 24 swg 0,5 mm
Suction cannulae
130 mm



WULLSTEIN SUCTION TUBE

29-960

Micro suction tube,

- 29-962** Suction tube with LUER hub
- 29-962** Micro suction cannula, 0.4 mm diam.
- 29-963** Micro suction cannula, 0.6 mm diam.
- 29-964** Micro suction cannula, 0.8 mm diam.



29-970

- 29-971** Suction handle, with finger cut-off, with LUER cone
- 29-972** 2 suction cannulae, 0.5 mm diam., with LUER cone
- 29-973** 2 suction cannulae, 0.7 mm diam., with LUER cone
- 29-974** 2 suction cannulae, 1.0 mm diam., with LUER cone
- 29-969** 2 suction cannulae, 1.3 mm diam., with LUER cone



ZOELLNER SUCTION TUBE

29-770
180 mm



BELLUCCI SUCTION TUBE

29-950

Suction set, complete, consisting of:

- 29-951** Suction tube with finger cut-off and stylet 155 mm



VERHOEVEN

diam.	diam.
29-763 16 swg	1,5 mm
29-764 18 swg	1,2 mm
29-765 20 swg	1,0 mm
29-766 22 swg	0,7 mm
29-767 24 swg	0,5 mm

Suction cannulae with LUER cone
100 mm



HOUSE

29-768

Suction adapter with finger cut-off,
LUER cone, LUER hub
55 mm



WULLSTEIN SUCTION TUBE

	diam.
29-239	1,5 mm
29-238	2,0 mm

Suction tube, cylindrical, silver, with stylet, olive with LUER hub
160 mm



MAGILL SUCTION TUBE - 245 MM

	diam.
29-906	3,0 mm
29-907	3,5 mm
29-908	4,0 mm



GILLIES SUCTION TUBE

	diam.
29-187	3,0 mm

with finger cut-off and stylet
170 mm



ROWBOTHAM SUCTION TUBE

29-672
with stylet
295 mm



FERGUSSON SUCTION TUBE

diam.	Length	Length	Length
	130 mm	160 mm	200 mm
1,5 mm	29-360	29-370	29-380
2,0 mm	29-361	29-371	29-381
2,5 mm	29-362	29-372	29-382
3,0 mm	29-363	29-373	29-383
4,0 mm	29-364	29-374	29-384
5,0 mm	29-365	29-375	29-385



Olive
diam. 6 - 9 mm

diam.	
29-360	2,0 mm
29-361	2,5 mm
29-362	3,0 mm
90° to curve, 250 mm	



olive diam. 6 mm

NEGUS SUCTION TUBE

29-177
Suction dissector
210 mm



olive diam. 6 - 10 mm

GORNEY SUCTION TUBE

29-168
Suction dissector
180 mm



wide

29-195 4 mm

29-196 5 mm

Suction raspatories, sharp, bayonet-shaped, with stylet
200 mm



LANDOLT SUCTION TUBE

wide

29-980 3,2 mm

Dissecting suction tube for hypophysectomie and ENT, sharp angled point
200 mm



Ø diam. 15 mm

FINSTERER SUCTION TUBE

29-867

320 mm



YASARGIL SUCTION TUBE

29-905
LUER-LOCK cone
200 mm



29-910
with finger cut-off and stylet



HOUSE SUCTION TUBE

	diam. suction tube	diam. irrigation tube
29-911	0,9 mm	0,7 mm
29-912	1,2 mm	0,9 mm
29-913	1,5 mm	1,2 mm

with finger cut-off and stylet, with LUER hub

